

HighTrack


High Pressure Low Compliant

Rapid Exchange PTCA Balloon Catheter

- ▼ Designed for use in highly calcified lesion and/or post dilation of stents
- ▼ Better option than Non-Compliant PET balloons as wall material is more resistant to puncture and burst, especially with highly calcified lesions.
- ▼ Highly Kink Resistant SPV Hypotube
- ▼ Low Compliant UHP (> 22 Atm) Balloon

Catheter Construction

Balloon Material	-	Low Compliant UHP Balloon Material
Balloon Compliance	-	Low Compliant
Proximal Shaft	-	SPV Hypotube
Distal Shaft	-	PEBAX
Max.Guidewire (Inch)	-	0.014"
Min.Guiding Catheter	-	5 Fr.
Marker Positioning	-	80 cm and 120 cm from distal tip
Diameters (mm)	-	2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0 mm
Lengths (cm)	-	10, 15, 20, 30 cm



Vasmed Technologies Limited

P O. Box 1356, Fujairah, United Arab Emirates
Ph.: +971 (0)9 2228862, Fax: +971 (0)9 2228863
Email: info@vasmed.ae

w w w . v a s m e d t e c h . c o m